

SIEMENS

MULTIMOBIL 2.5

SP

Disposal Instructions

© Siemens AG 2003

The reproduction, transmission or use of this document or its contents is not permitted without express written authority. Offenders will be liable for damages. All rights, including rights created by patent grant or registration of a utility model or design, are reserved.

English

Doc. Gen. Date: 11.03

Print No.: SPR8-X01.861.01.01.02

Replaces: n.a.

Multimobil 2.5

Med

Disposal Instructions
Environmental Protection Measures

Version: 01

Status: Released

Copyright © SIEMENS LTD.,MED INDIA

				Item Number	Version
Prepared by		Reviewed By		Released by:	
Name	Abhilash Kumar	Name	H.S.Usgaonkar	Name	Khandeparkar V.P
Dept.:	GW/Q	Dept.:	GW/R&D	Dept.:	GW / Q
Date:	22.01.2004	Date:	22.01.2004	Date:	28.01.2004
Signature:	Sd/-	Signature	Sd/-	Signature:	Sd/-

Index

1 HISTORY.....3

1.1 ALTERATION HISTORY.....3

1.2 REVIEW HISTORY3

2 GENERAL.....4

2.1 DISPOSAL LIST4

1 History

1.1 Alteration History

Version	Date	Prepared by / Dept.	Change and Reason for change
01	22-01-04	Abhilash Kumar	Not applicable being the first version

1.2 Review History

Version	Review date	Review protocol
---------	-------------	-----------------

2 General

Disposing of hazardous materials is as important when disposing of the system / equipments as it is when replacing individual items while the system is in use.

The list of hazardous materials in these instructions gives an overview of the components and assemblies requiring disposal. The disposal information and graphics will assist you in locating these hazardous materials on the system.

2.1 Disposal List

The items in the following table are to be disposed of by an authorised waste disposal company.

SI No.	Component Assembly	Location		Hazard/Hazardous Material	Qty.
1	X-ray tube assembly	Inside Single Tank	1/Fig1	Lead Oil X-ray tube implosion	1
2	Collimator	Below Single tank	2/Fig1	Lead	1
3	Printed Circuit Boards	Inside the control unit	3/Fig1	Epoxy resin on PC boards	1
4	Electrolytic Capacitors	Inside the control unit	4/Fig1	Electrolyte and electrodes	5




FIG. 1